

INVESTIGATION OF THE SEPARATION DEPENDENT FLUORESCENCE  
RESONANT ENERGY TRANSFER BETWEEN CDSE/ZNS QUANTUM DOTS BY  
NEAR-FIELD SCANNING OPTICAL MICROSCOPY

Pu Wang

Submitted to the faculty of the University Graduate School  
in partial fulfillment of the requirements  
for the degree  
Master of Science  
in the Program of Medical Biophysics and Biomolecular Imaging  
Indiana University

December 2009

Accepted by the Faculty of Indiana University, in partial  
fulfillment of the requirements for the degree of Master of Science.

---

Simon J. Atkinson, PhD, Chair

---

Ricardo S. Decca, PhD

Master's Thesis

Committee

---

Kenneth W. Dunn, PhD

---

Christoph A. Naumann, PhD

## ACKNOWLEDGMENTS

I would like to thank Dr. Ricardo Decca for all the help and support in this research.

I would like to thank Dr. Simon Atkinson for the help when I was in bad situation and support during the past two years.

I wish to thank Dr. Farbod Shafiei for his assistance and suggestions of the experiment. Without his hand-by-hand teaching, I could not complete this work.

I would like to thank Dr. Kenneth Dunn and Dr. Christopher Naumann for the course offering and knowledge sharing during the past two years.

Especially, this dissertation would not be completed without the love and support of my wife, Jing Ding.

## TABLE OF CONTENTS

LIST OF TABLES .....	vi
LIST OF FIGURES .....	vii
LIST OF ABBREVIATIONS .....	xii
ABSTRACT .....	xiii
1 INTRODUCTION .....	1
1.1 Quantum Dots .....	1
1.2 Förster Resonant Energy Transfer (FRET) .....	4
1.3 NSOM and Separation Dependence of Resonant Energy Transfer between CdSe/ZnS Nanocrystals .....	5
2 EXPERIMENTAL SETUP .....	9
2.1 Main Design .....	9
2.2 Near-field Scanning Optical Microscopy .....	11
2.2.1 Overview .....	11
2.2.2 Fiber Tips .....	12
2.3 Feedback Loop System .....	14
2.3.1 Tuning Fork .....	15
2.3.2 Quality of Engagement .....	19
2.3.3 Calibration of the distance between probe and cover glass .....	22

2.4	Poly(methyl methacrylate) (PMMA) and QDs Sample Preparation .....	23
2.5	CdSe/ZnS Quantum Dots for Resonant Energy Transfer .....	25
2.6	Spectrometer Setup .....	28
2.7	Photoluminescence Measurement.....	28
3	RESULTS AND DISCUSSION .....	32
3.1	Single Quantum Dot Localization by Coarse Scan.....	32
3.2	Small QDs-Large QDs Alignment by Fine Scan.....	36
3.3	Resonant Energy Transfer between CdSe/ZnS Nanocrystal Quantum Dots as Function of Separation .....	42
	APPENDICES .....	48
	REFERENCES .....	49
	CURRICULUM VITAE	

## LIST OF TABLES

Table

2.1	Table of calculated energy levels for the large QDs.....	26
-----	--	----

## LIST OF FIGURE

### Figure

1.1	Size dependent Photoluminescence (PL) spectra of CdSe QDs (from Benoit Dubertret and Hideki Ooba).....	2
1.2	Core/Shell structure of CdSe/ZnS Quantum Dots (QDs) .....	3
1.3	Schematic of NSOM probe and SIL (not to scale) with two resonant group of QDs deposited on probe and SIL. Red spheres representing large QDs and orange spheres representing small QDs. Gray hatched area is Al coating on the NSOM probe. The coated area prevents the excitation of dots which are not on the clear apex region.....	6
1.4	a) PL signal of small and large QDs as a function of separation. Horizontal, lateral and vertical axes show wavelength, separation between two groups of QDs and PL signal respectively. b) Quenching of PL signal of small QDs. Horizontal and vertical axis shows distance between two groups of QDs and the PL signal. The small QDs PL signal has been normalized by large QDs PL signal.....	7
1.5	Schematic of NSOM probe and cover glass (not to scale) with two resonant group of QDs deposited on probe and cover glass. ....	8
2.1	Schematic of experimental setup including near-field scanning optical microscope (NSOM) and spectrometer.. ....	10

2.2	a) SEM of 5-10nm Au coated fiber probe, magnification=35.000X. The arrow shows the flat apex. b) SEM of 15nm Cr plus 150nm Al coated fiber probe, magnification=150,000X .....	13
2.3	Schematic of feedback loop system. Direction of probe to cover glass approach and direction of probe/tuning fork oscillation are identified by black arrows.....	16
2.4	Schematic of RLC circuit representing tuning fork.....	17
2.5	Data and fitting of the real part of the current on a tuning fork as a function of frequency. The resonance frequency of the tuning fork is 32809Hz.....	18
2.6	a) Resonance curve of a tuning fork. b) Resonance curve of tuning fork with attached coated fiber probe.. .....	19
2.7	Response of the NSOM probe to a square wave set point. Vertical axis shows the amplitude of set-point (in red) as square wave and the NSOM probe amplitude (in black). Horizontal axis is time.....	20
2.8	Topography image of the micro beads. The right up corner shows the 2-D image and the left down corner shows the surface plot.....	21



2.9	Amplitude of NSOM probe as function of separation between NSOM probe and cover glass. Y-axis shows the resonance amplitude of the NSOM probe (measured by lock-in amplifier connected to the tuning fork and this value is the same as the setpoint). X-axis shows the separation between NSOM probe and cover glass.....	23
2.10	a) AFM image of the topography of the PMMA thin film. b) Line cut along dotted line in a). .....	24
2.11	Thickness of 4% PMMA on anisole thin film as function of the spinning rate. The thickness does not decrease if the speed is above 6000rpm. ....	25
2.12	Initial and final states of combination of small and large QDs. At initial state, the small QD exciton is in its ground state and there is no exciton in the large QD. In final state, there is not an exciton in the small QD and large QD exciton is in its third excited state.....	27
2.13	Side view of parts in the cryostat and out of the cryostat. The parts outside the cryostat have been designed to be used with objective for collecting the light and detecting it toward the spectrometer: 1-X-Y axis stage, 2-Z axis stage, 3-Mirror mount and mirror, 4-Objective tube, 5-Bracket for objective .....	30

2.14	sample PL spectra of QDs both from the probe and the cover slide. The green line shows the multi-peak fitting by Gaussian functions. These spectra are taken in 30s exposure time at 0.2mm entrance slit of spectrometer with notch filter applied in near field condition .....	31
3.1	Schematic shows FRET level measurement process. Area 1 is the intensity of the small QDs signal and area 2 is the intensity of the large QDs signal. The FRET level is calculated as area 1/area 2, which is the small QDs signal over the large QDs signal. Lower values indicate larger FRET interaction. ....	32
3.2	The upper row shows the 4x7 image of a) large QDs signal, b) small QDs signal and c) large QDs signal over small QDs signal. The lower row shows the pixel matrix corresponded with the images on the upper row. ....	34
3.3	a) number of pixels as a function of number of QDs we estimate. b) linear relationship between the pixel value and number of QDs. c) linear relationship between the pixel value and number of QDs with the un-normalized large QDs signal. The fitting parameters are also included. ....	35
3.4	PL spectra of small and large QDs showing one large QD signal.....	36
3.5	a) 7x7 fine scan image of small QDs signal (normalized by laser intensity. b) intensity value corresponding to the image. ....	37

3.6	a) $5 \times 7$ fine scan image of small QDs signal (normalized by laser intensity) scanning a blank cover glass. b) pixel value correspond to the image .....	39
3.7	a) $7 \times 7$ fine scan image of large QDs signal (normalized by laser intensity). b) Intensity corresponding to the image.....	40
3.8	a) The histogram of pixels value can represent three different populations of large QDs on the PMMA. b) fitting for the average pixel value which correspond to the three different population of large QDs .....	41
3.9	Schematic shows the possible situation happened in fine scan .....	42
3.10	a) PL spectra of small and large QDs as function of separation. Horizontal, lateral and vertical axes show wavelength, separation between the small QDs and large QDs. b) PL signal of small QDs normalized by the PL signal of large QDs. The horizontal axis shows the separation between small QDs and large QDs. The vertical axis shows the area of the PL spectra of the small QDs over the area of the PL spectra of the large QDs. c) small QDs signal when the probe approach to a blank cover glass.....	44
3.11	Transition probability as a function of the separation. The solid line shows the theoretical calculation by the former colleague. ....	45

## LIST OF ABBREVIATIONS

CCD: Charged Coupled Device

QDs: Quantum Dots

NSOM: Near-field Scanning Optical Microscope

FRET: Förster Resonant Energy Transfer

PL: Photoluminescence

## ABSTRACT

Pu Wang

### INVESTIGATION OF THE SEPARATION DEPENDENT FLUORESCENCE RESONANT ENERGY TRANSFER BETWEEN CDSE/ZNS QUANTUM DOTS BY NEAR-FIELD SCANNING OPTICAL MICROSCOPY

A Near-field Scanning Optical Microscope (NSOM) is used to study the resonant energy transfer between different size CdSe/ZnS quantum dots (QDs). The NSOM system is used to bring the small QDs which are 6 nm in diameter close to 8 nm diameter QDs which are embed with PMMA on a cover glass. The PMMA is used to prevent the 8 nm QDs from aggregation, which allows us to locate one dot on the cover slide and have the potential to get the interaction of two individual dots. A systematic methodology is used to localize a single QD on the cover glass and align the small and large QDs. Since the ground energy state of the small QDs match the excitation energy level of the large QDs. When the small dots get excited, part of the energy transfers to the large QDs. As the separation between small and large QDs is changed in near-field range (20-50nm), the transition probability is observed, indicating that the FRET level changes as a function of separation between small and large QDs. Possible future improvements are also discussed.

Simon J. Atkinson, PhD, Chair